
User manual

Rex

etac[®]
Creating Possibilities

78570 - 2023-11-01 - Version 5

sv en da de it nl no bg cs el es et fi fr hr is lt ms pl pt sk ar


www.etac.com


MD CE


A


8170 5010

8170 5020


8170 5040





1

8170 5030

8170 5050

B

| | | | | |
|---|---|---|---|---|
| 1 | 2 | 3 | 4 | 5 |
| 6 | 7 | 8 | | |

C

6

UDI Bar code GS1-128
(01)XXXXXXXXXXXXXXXX(11)YYMMDD(21)XXXXX

1 XXXXXXXX
XXXXXXX

2 **REF** XXXXXXXXXXXX
SN XXXXXX

CE MD Etac Supply Center AB
Långgatan 12
SE-334 33 Anderstorp
YYYY-MM-DD

3
4
5

D

| | | | | | | | | |
|---------------------------------------|---------------------------------------|------------------|------------------|------------------|--------------------------|---------------------------|---|---------------------------|
| | | | | | | | | |
| 600/700/850 mm (23.6"/27.6"/33.5") | 670/767/911 mm (26.4"/30.2"/35.9") | 155 mm (6.1") | 155 mm (6.1") | 130 mm (5.1") | 155x95 mm (6.1"x3.7") | 193x130 mm (7.6"x5.1") | 800 +70 -30 mm (31.5" +2.76" -1.18") | 52x55 mm (2.05"x2.17") |

| | |
|-----------------|----------------|
| | |
| *135 kg/297 lbs | 150 kg/330 lbs |

*85 cm=100 kg/220 lbs

| | | | | | | |
|-----------------|-------------------|-------------------|-------------------|-------------------|-------------------|--------------|
| | | | | | | |
| 8170 5010 | 8170 5020 | 8170 5040 | 8170 5030 | 8170 5050 | | |
| | | | | | | |
| 2.7 kg 6 lbs | 2.9 kg 6.4 lbs | 3.3 kg 7.3 lbs | 3,6 kg 7.9 lbs | 4,0 kg 8.8 lbs | Alu/ABS/PUR/PP/SS | PA6GF/TPE/SS |

E

1

800 mm

600 mm

2

| | |
|----|------------------|
| | B (mm) |
| 60 | 763 |
| 70 | 736 |
| 85 | 696 |

B

3

4

etac

5

etac

6

etac

7

Silicone

8

9

Silicone



10

11

Click!

Click!

12

F

1

2


3


G


1

2

3

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General information

Thank you for choosing an Etac product. In order to avoid damage during assembly, handling and use, it is important to read this manual and save it for future reference. You can also find it at www.etac.com. You can select your language via the "International" and "Local websites" link. Here you will also find other product documentation, such as prescriber information, pre purchase guide, and reconditioning instructions.

In the manual, the user is the person sitting in the chair. The carer is the person helping the user.

Device description

Rex is a foldable, wall-mounted toilet arm support, available in three different lengths.

Intended purpose

Rex (hereafter also referred to as "the device" or "the product"), is a medical device intended for alleviation of or compensation for a functional impairment due to an injury or disability. The device is designed to facilitate access to the toilet for people with restricted functional ability. The raised seat and/or arm supports intend to provide stability and support when sitting down and standing up before and after performing toilet needs.

Intended user group

The user group for the device is based on the individual's functional ability and not on a specific diagnosis, health condition or age. It is intended for individuals with a height of 146 cm or more, or a mass of 40 kg or more. Secondary users of the device are caregivers providing assistance and clinicians/technicians who set up the device.

Intended environment

The device is intended for indoor use at home environments or institutions and is suitable for use in bathrooms, but not swimming pools or similar corrosive environments.

Intended application

The device is intended for short-term and long-term use and can be applied several times per day. The device is intended to be used in contact with intact skin. The device is intended for refurbishment and reuse.

Expected service life

The expected service life is 10 years. For complete information regarding the service life of the device, see www.etac.com.

Indications

Indication for use is disability, including, but not limited to, an individual with fluctuating balance issues, physiological and/or functional impairment, and/or a history of falls.

Contraindications

There are no known contraindications.

Warnings

Warnings describing a risk element for a specific action or setting of the device can be found in the relevant section.

Mounting on walls with insufficient load-bearing capacity or with unsuitable fixtures can affect user safety. Mounting of the device and assessment of the wall's load-bearing capacity must be carried out by a qualified person.

A defective device must not be used.

Ignoring installation or assembly instructions may result in injury.

Do not stand on the device

Declaration of conformity

The device meets the requirements of the Medical Device Regulation (EU) 2017/745.

The device is tested and meets the requirements of the EN ISO 10993-1, SS EN 12182.

Pre-purchase information

Pre-purchase information can be found at www.etac.com.

Adverse event

In case of an adverse event occurred in relation to use of the device, it should be reported to your local dealer and the national competent authority in a timely manner. The local dealer will forward information to the manufacturer.

Customization

is everything that goes beyond the manual's instructions and settings. Device specially adapted by the customer must not retain Etac's CE marking. Etac's warranty expires.

If there is any uncertainty regarding the execution of the customization, consult Etac.

Device coupled with another device may not retain either device's CE mark unless combination agreements exist. Contact Etac for current information.

Device description [Figure A](#)

1. Content of the delivery **2.** Supporting legs **3.** Wall bracket **4.** Arm support **5.** Adjustable ferrule **6.** Product label **7.** Item no label

Symbols [Figure B](#)

Symbols in the manual and on the device:

1. Warning, precaution or limitation. **2.** Useful advice and tips. **3.** Material for recycling. **4.** Maximum user weight (see technical data) **8-9.** Cleaning (see maintenance) **7.** Read instructions for use **8.** Product weight

Explanation item no label: [Figure C](#)

1. Product name **2.** Product description **3.** Serial number **4.** Item number **5.** Date of manufacture **6.** Barcode acc. to GS1-128 GTIN-14 and serial number*

*The device's date of manufacture can be read from the barcode on the device. The number 11 is shown below the barcode in brackets. The number combination after these brackets is the date of manufacture.

Technical data [Figure D](#)

Installation options

May be mounted on walls with sufficient load-bearing capacity only. Please ensure that suitable fixtures are used. For mounting in concrete, plywood or wooden joists, we recommend the screw set accessory.

Installation [Figure E](#)

Use [Figure F](#)

Maintenance

Clean the device with a solvent-free detergent

which has a pH of 5-9, or with a 70% disinfectant solution. Can be decontaminated at max. 85 °C if so required.

Constituent materials are corrosion-resistant. 5-year guarantee against material and manufacturing defects. For terms and conditions, see www.etac.com.

Storage

The device should be stored indoors in a dry place at a temperature above 5 °C. If the device has been stored for a long time (more than four months), its function must be checked by an expert before use.

Accessories [Figure G](#)

1. Toilet paper holder **2.** Storage basket **3.** Screw kit

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Etac Ltd.
Promefa AG