



Etac Optima toilet arm supports

Designed to assist users and carers with every day routines.

 **etac**[®]
Creating Possibilities

Etac OptimaL toilet arm supports

Modern and elegantly designed arm supports.

Function and simplicity in an elegant combination

OptimaL is an elegant wall mounted support rail. The leg secures the stability and the small wall attachment enables mounting close to the toilet. OptimaL stays safely folded until it is flipped down with the grey handle.

Accessories

See our wide range of accessories at the end of this brochure



compact folding

Description	Item no.	
Optima L	80303006	✓

Dimensions

Wall to arm supports	73.5 cm 29"
Floor to arm supports	80 cm 31½"
Height, floor to top of arm supports when folded up	149 cm 58¾"
Wall bracket	10.4 x 11.0 cm 4" x 4¼"

Material

Wall bracket, arm supports tubing, supporting leg	Powder-painted steel tubing
Details	Polypropylene

Weight

3.1 kg 6.8lbs



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.

Etac Optima toilet arm supports

Optima toilet arm supports provides a firm and comfortable grip and folds smoothly to the wall.

Comfortable and well-fitting

Optima arm support cover has a softly rounded end which provides a comfortable grip. The supporting leg has an adjustable foot which ensures a perfect floor contact. The large wall plate gives a firm fixation to the wall.



Description	Item no.	
Optima	80303001	✓

Dimensions

Wall to arm supports:	71.5 cm 28¼"
Height arm supports:	80 cm 31½"
Height folded arm supports (from floor):	147.5 cm 58"
Wall bracket:	16.5 x 30.4 cm 6½ x 12"

Material

Wall bracket:	Powder-painted aluminium
Arm supports tubing, supporting leg:	powder-painted steel tubing
Arm supports washer:	Polypropylene

Weight

2.4 kg 5.3 lbs



max 85°C / max 185°F
for maximum of 3 minutes

max 130 kg / max 287 lbs



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.



compact folding

Etac Optima Lockable toilet arm supports

Optima lockable offers extra security and support.

Unique locking mechanism provides extra stability

Optima lockable locks both in the folded down and up position. In the locked folded up position, it provides the user support in a standing position. The unique locking mechanism is easy to operate with one hand. The arm supports cover has a softly rounded end which provides a comfortable and secure grip.

Accessories

See our wide range of accessories over



Optima Lockable is available in two different lengths.

Description	Weight	Item no.
Optima Lockable	2.4 kg 5.3lbs	80303004

Dimensions

Length arm support	71.5 cm 28¼"
Wall bracket	16.5 x 30.4 cm 6½" x 12"

Material

Wall bracket	Powder-painted aluminium
Arm supports tubing	Powder-painted steel tubing
Locking device	Polyamide
Arm supports washer	Polypropylene



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.



compact folding

Etac Optima accessories



Description Attachment for signal device
Item no. 80303027

Suitable for
Optima
Optima Lockable



Description Floor pillar
Item no. 80303028

For floor-mounting, positions arm supports height at 80 cm 31½"

Suitable for
Optima
Optima Lockable



Description Supporting leg
Item no. 80303040

For relieving strain on walls.

Suitable for
Optima Lockable only



Description Extension, 10 cm 4"
Extension, 20 cm 7¾"
Extension, 30 cm 11¾"
Item no. 80303030
80303031
80303032

Suitable for
Optima
Optima Lockable



Description Toilet paper holder
Item no. 83030001
Complete, fitted to the arm supports.


Suitable for
OptimaL
Optima
Optima Lockable





 [®]
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item