

# Etac Turner Pro

Safe transfers that are kind to a carers back



## **Etac Turner PRO**

Patient turning with standing support

## The importance of safe transfers

People in need of assistance often feel unsafe in transfer situations, and caregivers are sub-ject to strain when assisting with the transfer. The risk of accidents and insecurity from both caregiver and the assisted person may result in fewer transfers.

The Turner PRO allows both the assisted person and the caregiver to feel safe and the transfer performed with ease.

## Easy to use

Turner PRO is the perfect tool for hospitals, care homes and in the home. The intiutive design makes it easy to use, the only adjust-ments needed are handle height and height of leg support.

The Turner PRO works well on different flooring materials including carpets.



A safe transfer which allows eye contact and good communication.

## **Turner PRO**

Code 16090106

#### **Dimensions**

Handle height: 78-120 cm / 30.7"-49.2" Footplate diameter: 40 cm / 15.7"Footplate height: 20 mm / 0.8"Lower leg support width: 38 cm / 15"Lower leg support height: 15 cm / 5.9"

#### Material

Cast aluminum, powder coated steel, chrome 3

## Compliant with

MDD 93/42/EEC EN 12182

## Total weight

8 kg / 17.6 lbs

## Weight

Max user weight 200 kg / 441 lbs







Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decomtaminated at max  $85^{\circ}$ C /  $185^{\circ}$ F is so required.

## **Features**



Compact

Easy to transport. Easy to grip. Light, evenly balanced weight.



## **Angled leg support**

Soft and comfortable, provides lateral support. Easy to adjust to whatever height is needed.



### Contrasting coloured handle

The contrasting colour give the handle extra visibility. Offers many grip possibillities for both carer and user.



### **Functional base plate**

Cast aluminum base plate with integrated anti slip function. Very easy to turn. The wheels make Turner PRO easy to transport indoors.



Design, function and safety are extremely important to us, and this is presented here in the form of a sit to stand aid. The scandinavian design gives you maximum comfort and a more stylish environment.



Turner PRO has a small footprint which makes it the perfect choice when transferring in a tight space, for example between a wheelchair and toilet.



Encourages natural standing and sitting movement for the user. Enables a good posture and minimises the risk of injury for the carer when transferring.



## Intended Use

Turntable with support handle

- Transfer from bed to wheelchair / shower chair
- · Transfer from wheelchair to toilet
- Transfer from chair to wheelchair / shower chair

Go to **novis.com.au/turner-pro** or scan the OR code on the back cover to find out more.



"Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances"









