









# Multiple Positioning and Reduced Shear

The BetterLiving Community Bed is incredibly practical yet versatile bed to suit a range of home care needs. Operates as a four part, profiling bed with electric high-low height adjustment, back rest, knee break and Trendelenburg tilt – to deliver multiple positioning options to suit every need.

The Auto Contour feature simultaneously raises both the head and knees, which combined with auto-regressive backrest helps prevent painful skin shearing.

# Quality Design

Powered by a high quality Linak electric drive and constructed from durable powder coated steel. The frame includes a 200mm Length Extension and removable mattress retainers.

Designed with built-in storage features for the Length Extension and Transport Bracket which minimises misplacement when not in use.

Fitted with four high quality individually lockable castors for easy maneuverability.

### Easy, Safe Operation

Lightweight, water resistant hand controller with soft touch buttons for raising, lowering and positioning, including one touch auto contouring.

Convenient inbuilt hanging hook for easy storage and access. Features lockout function, ideal for patients incapable of operating the bed due to cognitive impairment.

Each function can be individually locked out using the key, which is required to unlock and reactivate functionality.



# Community Bed

SIZE		Single	King Single	Bea
CODE		BE1005	BE2005	
DIMENSIONS	LENGTH	2170 mm to 2370 mm (Extended)	←	
	WIDTH	995 mm	1220 mm	
BED FRAME WEIGHT		74 kg	84 kg	
MATTRESS SIZE RECOMMENDED	WIDTH	880 mm	1050 mm	
	LENGTH	1980 mm	←	
	HEIGHT	170 maximum	←	
SAFE WORKING LOAD		220 kg	←	
MAXIMUM PATIENT WEIGHT		180 kg	←	
WARRANTY		5 year bed frame / 2 year electrical / 1 year castor		
FRAME MATERIAL		Powder Coated Steel	←	
HEIGHT RANGE		330 mm to 730 mm	←	
BACKREST ANGLE		0° to 70°	←	
KNEE BREAK ANGLE		0° to 32° (4 Section Bed)	←	
TRENDELENBURG / REVERSE TILT		0° to 10°	←	
COMPLIANCES		IEC 60601-1:2005+A1:2012 [Medical Electrical Equipment Part 1]		
		IEC 60601-2-52: 2012 [Medical Beds Safety and Essential Performance]		
		IEC 61000-3: 2008 [EMC Electromagnetic compatibility]		
		IEC 61000-4: 2009 [EMC Electromagnetic compatibility]		
		CISPR 11: 2004+A2: 2006, Group 1 Class B [EMC Electromagnetic compatibility]		

#### Accessories



BE6003 Side Rail, Three Quarter Length (Sold as Pair)



BE4017 Transport Bracket (Sold as Pair)



BE8881 BE8882
Single King Single
Slip Cover, Grey
(Sold as a pair for Head and Foot)



BE9002 Grab Bar



transport bracket BE4017



For more information: www.novis.com.au/ product/betterlivingcommunity-bed