

# Tapered Edge Non-Slip Bedside Safety Mat



Even though the role of bedrails in falls prevention is well established, there may be circumstances where the use of such equipment may not be appropriate, such as:

- For a patient who could be independently mobile without bedrails
- Patients who are distressed by the use of bedrails
- Severely confused patients who may attempt to climb over the bedrails

In such patients, who are at risk of falls from beds, alternative equipment, such as low/floor beds and Safety Mats, are indicated.

FEATURE	FUNCTION	BENEFIT
<b>Tapered Edge</b>	Reduced trip hazard	Minimises falls risk
Nylon Material	Low friction	Easy to slide under bed without the need to bend
Textured Base	Non-slip when pressure applied	Will stay in position if stood on
<b>Reflective piping</b>	Increased visibility at night	Reduces the risk of tripping
<b>Three hinged sections</b>	Foldable	Easy daytime storage
<b>Instruction Tag</b>	Outlines product care requirements	Assists with improved usage and care of product



# Tapered Edge Non-Slip Bedside Safety Mat



This Tapered Edge Safety Mat can be placed on either side of the bed as a protective cushion to help prevent injury for patients at risk of rolling out of bed. Especially useful where side rails cannot be used.

The non-slip feature helps to minimise the risk of injury should the mat be inadvertently stood on by patient or carer. The unique textured underside prevents the mat slipping when downward pressure is applied however still allows the mat to be repositioned with the foot for easy manoeuvrability.

The Safety Mat includes a reflective strip to minimise trip hazards in low light visibility. Comprises three sections of safety foam covered by hinged nylon for easy folding and storage.

<b>CODE</b>	BL7274
<b>THICKNESS</b>	50 mm
<b>WIDTH</b>	750 mm
<b>LENGTH</b>	1750 mm
<b>CONFORMS TO</b>	ARTG 155475
<b>WARRANTY</b>	2 year