

Etac My-Loo

Raised toilet seat, with brackets

(oetac)

Etac My-Loo

Raised toilet seat, with brackets

Together with experienced industrial designers, we have developed the new raised toilet seat with brackets – Etac My-Loo. A product that gives concepts such as hygienic, safe and appealing a new meaning. With its gently sweeping lines and excellent ergonomics, it offers extra-comfort and a design ideally suited to the modern bathroom.

Designed for function and comfort

Etac My-Loo has been developed to blend into the bathroom environment and has a design that gives optimal comfort and accessibility. The gentle lines and the large opening make My-Loo suitable for most users. The seat raiser is offered in two heights, with or without a lid.

Easy to keep clean and fresh

My-Loo raised toilet seat, with brackets, consists of a small number of parts made of durable materials. It is easy to wipe off with a soft cloth. Should separate cleaning be required, it is easy to remove the seat and equally simple to put it back in place.

Safe, quick installation

It takes only a few seconds to install Etac My-Loo! You fix it in place using two specially designed brackets which are easily adapted to different toilet shapes. The extra-large friction surfaces and flexible material of the brackets provide a firm attachment to the porcelain.

Put the raised toilet seat on the toilet and push the brackets into place. Click, and the installation is finished!



Etac My-Loo raised toilet sea with brackets	rt, Item no.	Design
with bidokets		DESIGN
6 cm	80301520	AWARD DOLLES
10 cm	80301521	red <mark>dot</mark> award 2019
6 cm with lid	80301522	
10 cm with lid	80301523	Myra Industriell Design AB
Weight		Maintenance
6 cm: 0.8 kg	10 cm: 1.1 kg	
6 cm with lid: 1.2 kg	10 cm with lid: 1.4 kg	Alk pH 70% 5-9 23 min MAX 190 kg 418 lbs
Maximum user weight		
190 kg		Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH
100 Kg		between 5 and 9, or with a 70% disinfectant solution. The product can be safely cleaned
Material		in a cabinet washer-disinfector at 85°C for 3 minutes.
Polypropylene and TPE		

2 Etac My-Loo Etac | Creating Possibilities

Features



Easy to keep clean and fresh

The clean, smooth surfaces are easy to clean with a soft cloth. Should separate cleaning be required, the seat is easy to remove and equally simple to put back in place.



Designed for function and comfort

The soft shape of the seat gives safe support and high comfort. An extra large opening and generous recesses facilitate accessibility during personal hygiene.



Safe, quick installation

The easily attached brackets with wide adjustability and large friction surfaces make My-Loo at all times safely fixed to, and easily mounted on, any toilet porcelain. Entirely without screws.



Variants

My-Loo raised toilet seat, with brackets, is available in two different heights, 6 and 10 cm, with or without lid.



Dimensions

Outer width x depth 36 cm x 40 cm
Opening width x depth 23 cm x 36 cm
Seat height 6 cm, 10 cm
Distance between brackets 29 cm - 38,5 cm













3

311 X 40 CH1 23 CH1 X 30



Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Our heart lies in the solutions that optimise quality of life for the individual, their family and caregivers.

For the latest news and continuously updated product information – visit: www.etac.com

Etac AB
Kista Science Tower
SE-164 51 Kista, Sweden
Tel +46 371 - 58 73 30
info@etac.se www.etac.com

