

# Compatibility Matrix - 2701 Falls Monitor with SafeTime

Suitable for private residences, hospitals and care homes, its versatility turns a basic room into a highly optioned care and monitoring facility offering multiple solutions in one simple to operate unit.



- Durable ABS plastic enclosure with soft touch removable protection cover and Caregiver Reset Key
- Amber LED visual alarm indicator
- Caregiver secure reset option
- Battery operated or optional external plug pack
- 6 - 12 months battery life (nominal)
- Easy to access battery compartment for staff
- Easily connected to a nurse call system, selectable alarm tones, 0, 2 sec, and 15min selectable alarm delay settings
- Easy thumbwheel adjustable alarm volume
- "Pause" feature allows caregiver to attend to the patient without initiating an alarm (will auto reset after 3 minutes).
- Alerts for lost nurse call connection, lost AC power & low battery
- Change Pad alarm feature when SafeTime pads used






### What's in the box:

- 1 x 2701 Standard Falls Monitor
- 1 x Quick Start Guide
- 1 x 2642 Caregiver Reset Key
- 1 x Durable Soft Protective Cover
- 3 x Alkaline AA Batteries

|   | Code | Product                               | Self Reset | Description   |
|---|------|---------------------------------------|------------|---|
|  | 3516 | Premium Under Mattress Bed Pad        | Y          | Positioned UNDER mattress. The monitor alarms when getting out of bed. The monitor auto resets when back in bed.            |
|  | 2626 | Premium Bed Pad                       | Y          | Positioned ABOVE mattress. The monitor alarms when getting out of bed. The monitor auto resets when back in bed.            |
|  | 2628 | Premium Wheelchair / Chair Pad        | Y          | Positioned on chair / wheelchair. The monitor alarms when getting out of chair. The monitor auto resets when back in chair. |
|  | 2633 | Premium Heavy-Duty Non-Slip Floor Mat |            | Positioned strategically in an area of risk. The monitor alarms when stood on.  |

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

|   | Code | Product                                   |   | Description  |
|---|------|---|---|--|
|    | 2623 | SafeTime Bed Pad - 1 Year                 | Y | Positioned above mattress. The monitor alarms when getting out of bed. The monitor auto resets when back in bed. The monitor alarms when pad expires.              |
|    | 2622 | SafeTime Bed Pad - 90 days                | Y | Positioned on a chair / wheelchair. The monitor alarms when getting out of bed. The monitor auto resets when back in bed. The monitor alarms when pad expires.     |
|    | 2621 | SafeTime Wheelchair / Chair Pad - 1 Year  | Y | Positioned on a chair / wheelchair. The monitor alarms when getting out of chair. The monitor auto resets when back on chair. The monitor alarms when pad expires. |
|   | 2620 | SafeTime Wheelchair / Chair Pad - 90 days | Y | Positioned on a chair / wheelchair. The monitor alarms when getting out of chair. The monitor auto resets when back on chair. The monitor alarms when pad expires. |
|  | 2742 | Spare Soft Cover for Monitors             |   | Soft touch removable protection case provides extra protection for the 2701 Monitor.   |
|  | 2406 | cura1 6Vdc Plug Pack                      |   | Optional Plug Pack.  |
|  | 2642 | Caregiver Reset Key                       |   | Resets alarms on Falls Monitor.  |
|  | 3504 | ActiveCare Mounting Bracket               |   | Monitor clips into Mounting Bracket. Mounting accessories included.  |
|  | 2636 | Easy Release Wheelchair / Chair Seatbelt  |   | Activates an alarm if buckle is released.  |

|   | Code | Product                                    |   | Description   |
|---|------|--|---|---|
|    | 2643 | Nurse Call Cable 6.35 Mono to 6.35 Mono 2m | Y | Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.) |
|    | 2646 | Nurse Call Cable 6.35 Mono to 6.35 Mono 3m |   | Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.) |
|    | 2493 | Nurse Call Cable 6.35 Mono to 3.5 Mono 2m  |   | Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.) |
|   | 2461 | Nurse Call Cable 6.35 Mono to 6.35 T&S 2m  |   | Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.) |
|  | 2460 | Nurse Call Cable 6.35 T&S to 6.35 T&S 2m   |   | Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.) |

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.