

Compatibility Matrix – 2544 Multiport Expander






Suitable for nursing homes and hospitals where nurse call systems are installed.

- Works with up to 5 wireless devices
- Operating range 8 – 12 metres
- 6.35mm N/O plug for connection to N/C system
- 6.35mm socket for connecting a nurse call pendant
- RJ10 socket for connecting hardwired pads & mats
- 3 x AAA Batteries
- 2406M 6V Plug Pack (optional)
- Easy to program
- Low Battery and Pad Lost indication
- Available in Tip and Ring configuration – 2544TR
- Available in Ring and Sleeve configuration – 2544RS







What's in the box:

- 1 x 2544 Multiport Expander
- 1 x Quick Start Guide
- 3 x High quality AAA batteries

	Code	Product	Self Reset	Description
	3515	Cordless Under Mattress Bed Pad	Y	Located UNDER MATTRESS. Sends OUT OF BED message when getting out of bed, triggering the nurse call system.
	2629	Cordless Bed Pad	Y	Located on TOP OF THE MATTRESS. Sends OUT OF BED message when getting out of bed, triggering nurse call system.
	2630	Cordless Chair Pad	Y	Sends OUT OF CHAIR message when getting out of chair, triggering the nurse call system.
	2624	Cordless Floor Mat		Sends message when stepped on to Multiport Expander triggering nurse call system.
	2719	Cordless Pad Transmitter		Required by Cordless Pads/Mats, sends a signal to the 2700.

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

	Code	Product		Description
	2626	Premium Bed Pad	Y	Sends OUT OF BED message through RJ10 Cord connected to Multiport Expander, triggering the nurse call system.
	3516	Premium Under Mattress Pad	Y	Sends OUT OF BED message through RJ10 Cord connected to Multiport Expander, triggering the nurse call system.
	2628	Premium Chair Pad		Sends OUT OF BED message through RJ10 Cord connected to Multiport Expander, triggering the nurse call system.
	2633	Corded Floor Mat Black/Grey		Sends OUT OF BED message through RJ10 Cord connected to Multiport Expander when person stands on floor mat, triggering the nurse call system.
	2588	Personal Emergency Transmitter (PET)		Sends message when button is pressed, triggering the nurse call system.
	2721	iBeam		Sends alarm message when triggered. Features a narrow detection zone to reduce false alarms. Triggers when a patient breaks the beam.
	2638	Nurse Call Pendant		Simple operation. Plugs into the 2544 Multiport expander. N/O push button. 2 metre flexible cord.
	2406M	cura1 6V Plug Pack		Plug Pack with 1.3mm connection tip.

	Code	Product	Description
	2558	Wireless Waterproof Call Button	Sends message when ORANGE button is pressed, triggering nurse call system.
	3502	ActiveCare Wireless Handset	Sends message when RED button pressed, triggering nurse call system.
	3503*	Simple Push Button Call Point*	Sends message when button pressed, reset at console.
	3503	ActiveCare Magnetic Reed Switch	Sends message to Multiport Expander when button is pressed or door is opened, triggering nurse call system.

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.