

# Compatibility Matrix – 2549 Cordless Multiport Expander with Bracket Mount






Suitable for nursing homes and private public hospitals where nurse call systems are installed.

- ☐ Works with up to 5 wireless devices
- ☐ Operating range 8 – 12 metres
- ☐ RJ12 socket for connecting to N/C system
- ☐ 3 x AAA Batteries
- ☐ Optional 2406M 6V Plug Pack
- ☐ Easy to program
- ☐ Low Battery and Pad Lost indication














### What's in the box:

- 1 x 2549 Multiport Expander Bracket Mount
- 1 x 3504 ActiveCare Bracket
- 1 x Quick Start Guide
- 3 x AAA Batteries
- 1 x Lead to connect to N/C system



|   | Code | Product                              | Self Reset | Description  |
|---|------|--------------------------------------|------------|--|
|  | 3515 | Cordless Under Mattress Bed Pad      | Y          | Located UNDER MATTRESS. Sends OUT OF BED message when getting out of bed, triggering the nurse call system.        |
|  | 2629 | Cordless Bed Pad                     | Y          | Located on TOP OF THE MATTRESS Sends OUT OF BED message when getting out of bed, triggering the nurse call system. |
|  | 2630 | Cordless Chair Pad                   | Y          | Sends OUT OF CHAIR message when getting out of chair, triggering the nurse call system.                            |
|  | 2624 | Cordless Floor Mat                   |            | Sends message when stepped on to Multiport Expander, triggering the nurse call system.                             |
|  | 2588 | Personal Emergency Transmitter (PET) |            | Sends message when RED button pressed, triggering the nurse call system.   |

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

|   | Code  | Product   | Description  |
|---|-------|---|--|
|    | 2721  | iBeam   | Sends alarm message when triggered. Features a narrow detection zone to reduce false alarms. Triggers when a patient breaks the beam.  |
|    | 2723  | EzyStand  | EzyStand can be installed adjacent to a bed or chair, or across a doorway. With our sturdy, adjustable-height EzyStand, precise positioning of the iBeam could not be simpler. |
|    | 2595  | Wristband Alarm Transmitter with Man Down           | Initiates alarm when RED button pressed, detects patient falls and automatically sends alarm. Reset on Falls Monitor.  |
|    | 2406M | cura1 6V Plug Pack                                  | Plug Pack with 1.3mm connection tip.   |
|   | 2558  | Wireless Waterproof Call Button                     | Sends message when ORANGE button is pressed, triggering nurse call system.   |
|  | 3502  | ActiveCare Wireless Handset (Includes wall bracket) | Initiates alarm when RED button pressed. Reset on Falls Monitor.   |
|  | 3503* | Simple Push Button Call Point*                      | Initiates alarm when pressed. Reset on Falls Monitor.  |
|  | 3503  | ActiveCare Magnetic Reed Switch                     | Initiates alarm when door is opened. Reset on Falls Monitor.   |
|  | 4160  | Curly Cord RJ12 to 6.35 R/A Plug St 300mm T&S       | Connects the 2459 directly to a nurse call system.   |
|  | 4161  | Curly Cord RJ12 to 6.35 R/A Plug St 300mm T&R       | Connects the 2459 directly to a nurse call system.   |
|  | 4162  | Curly Cord RJ12 to 6.35 R/A Plug St 300mm R&S       | Connects the 2459 directly to a nurse call system.   |
|  | 4163  | Curly Cord RJ12 to 3.5 R/A Plug St 300mm T&S        | Connects the 2459 directly to a nurse call system.   |
|  | 4164  | Curly Cord RJ12 to 3.5 R/A Plug St 300mm T&R        | Connects the 2459 directly to a nurse call system.   |