

Compatibility Matrix - 2701 Falls Monitor with SafeTime

Suitable for private residences, hospitals and care homes, its versatility turns a basic room into a highly optioned care and monitoring facility offering multiple solutions in one simple to operate unit.

- Durable ABS plastic enclosure with soft touch removable protection cover and Caregiver Reset Key
- Amber LED visual alarm indicator
- Caregiver secure reset option
- Battery operated or optional external plug pack
- 6 – 12 months battery life (nominal)
- Easy to access battery compartment for staff
- Easily connected to a nurse call system, selectable alarm tones, 0, 2 sec, and 15min selectable alarm delay settings
- Easy thumbwheel adjustable alarm volume
- "Pause" feature allows caregiver to attend to the patient without initiating an alarm (will auto reset after 3 minutes).
- Alerts for lost nurse call connection, lost AC power & low battery
- Change Pad alarm feature when SafeTime pads used



What's in the box:

1 x 2701 Standard Falls Monitor

1 x Quick Start Guide










1 x 2642 Caregiver Reset Key






1 x Durable Soft Protective Cover

3 x Alkaline AA Batteries

	Code	Product	Self Reset	Description
	3516	Premium Under Mattress Bed Pad	Y	Positioned UNDER mattress. The monitor alarms when getting out of bed. The monitor auto resets when back in bed.
	2626	Premium Bed Pad	Y	Positioned ABOVE mattress. The monitor alarms when getting out of bed. The monitor auto resets when back in bed.
	2628	Premium Wheelchair / Chair Pad	Y	Positioned on chair / wheelchair. The monitor alarms when getting out of chair. The monitor auto resets when back in chair.
	2633	Premium Heavy-Duty Non-Slip Floor Mat		Positioned strategically in an area of risk. The monitor alarms when stood on.

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

	Code	Product		Description
	2623	SafeTime Bed Pad - 1 Year	Y	Positioned above mattress. The monitor alarms when getting out of bed. The monitor auto resets when back in bed. The monitor alarms when pad expires.
	2622	SafeTime Bed Pad - 90 days	Y	Positioned on a chair / wheelchair. The monitor alarms when getting out of bed. The monitor auto resets when back in bed. The monitor alarms when pad expires.
	2621	SafeTime Wheelchair / Chair Pad - 1 Year	Y	Positioned on a chair / wheelchair. The monitor alarms when getting out of chair. The monitor auto resets when back on chair. The monitor alarms when pad expires.
	2620	SafeTime Wheelchair / Chair Pad - 90 days	Y	Positioned on a chair / wheelchair. The monitor alarms when getting out of chair. The monitor auto resets when back on chair. The monitor alarms when pad expires.
	2742	Spare Soft Cover for Monitors		Soft touch removable protection case provides extra protection for the 2701 Monitor.
	2406	cura1 6Vdc Plug Pack		Optional Plug Pack.
	2642	Caregiver Reset Key		Resets alarms on Falls Monitor.
	3504	ActiveCare Mounting Bracket		Monitor clips into Mounting Bracket. Mounting accessories included.
	2636	Easy Release Wheelchair / Chair Seatbelt		Activates an alarm if buckle is released.

	Code	Product		Description
	2643	Nurse Call Cable 6.35 Mono to 6.35 Mono 2m	Y	Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.)
	2646	Nurse Call Cable 6.35 Mono to 6.35 Mono 3m		Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.)
	2493	Nurse Call Cable 6.35 Mono to 3.5 Mono 2m		Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.)
	2461	Nurse Call Cable 6.35 Mono to 6.35 T&S 2m		Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.)
	2460	Nurse Call Cable 6.35 T&S to 6.35 T&S 2m		Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.)

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.