

cura¹ Fall Prevention Products



We have a range of products to assist with High Falls Risk patients /residents.

They can either cushion the fall from low lying beds or detect movement before they leave the bed.



PreAlert Bed Pad




Sensored Floor Crash Mat



Solutions

See our range of solutions below to either cushion the fall from low lying beds or early movement detection sensors.

Please see overleaf for assistance with ordering or contact us for our pricelist or further help.

Type	Picture	Best Features	Monitors Under 40kg	Alarms Both Sides of The Bed	Hardwired Options	Wireless Options
Sensored Floor Crash Mat		Mats are used beside low lying beds. Will cushion fall, diminishes injuries. Easily folded for storage. Simple operation and maintenance.	Yes	No	Yes	No
Non - sensed Floor Crash Mat		Mats are used beside low lying beds. Will cushion fall, diminishes injuries. Easily folded for storage.	-	-	-	-
PreAlert Bed Pad		Bed Pad gives an early warning when a person sits or rolls to the edge of the bed. It has a delay zone in the centre area to minimise false alarms. Has a 15 – month warranty & inbuilt transmitter.	No	Yes	No	Yes







Novis Healthcare
www.novis.com.au
1300 738 885


HOW TO ORDER:

Choose the monitor & then the pads that best suits your patients / resident's requirements. Don't forget your transmitter for wireless and leads for hardwired options.


Wireless Option Monitors:

Picture	Code	Description	More information
	2544	Multiport Expander	Can be programmed with up to 2 pads & 3 other transmitters. Silent within room. No reset button, staff cancel at call point only. Low battery & pad lost notifications.
	2549	Multiport Expander Bracket Mount	As above. Monitor comes with included changeable curly cord nurse call lead. Can be used with different nurse call connections or hard to reach call points.
	2702	3 Channel Wireless Monitor	Can be programmed with up to 3 wireless pads. Standalone or through nurse call system (requires connection lead). Can be silent within room or emit an audible alarm for extra urgency.
	2700	6 Channel Wireless Monitor	Can be programmed with up to 6 wireless pads. Standalone or through nurse call system (requires connection lead). Can be silent within room or emit an audible alarm for extra urgency.



Wireless Programmable Transmitter:

Picture	Code	Description	Best Feature	Transmitter Required
	2670	Wireless PreAlert Bed Pad	Alerts when they move to the edge of the bed	-

Hardwired Option Monitor:




Picture	Code	Description	More information
	4017	Economy Falls Monitor	Monitor for Economy corded pads, can connect to one pad. Adjustable alarm tone & volume for extra urgency. Standalone or through nurse call system (lead required). Low battery indicator. Monitor emits an audible alarm for extra urgency.

Hardwired Pads

Picture	Code	Description	Best Feature	Break Cord
	4114	Sensored Crash Mat	Will cushion fall, diminishes injuries. Requires no monitor, lead connects direct to call point. Lead will need to be changed if using the monitor.	4150 Included
	4113	Non - Sensored Crash Mat	Will cushion fall, diminishes injuries.	-

Connection Leads/ Break Cord:

If using a 2702 3 Channel Monitor, 2700 6 Channel Monitor or 4017 Economy Monitor and you would like to connect this to a call point you will require a nurse call connection lead.

Picture	Code	Description	More information
	2461	Nurse Call Connection Lead	6.35mm mono to 6.35mm stereo (TS) – 2m
	2643	Nurse Call Connection Lead – suits most systems	6.35mm mono to 6.35mm mono (TS) – 2m
	2646	Nurse Call Connection Lead – suits most systems	6.35mm mono to 6.35mm mono (TS) – 3m
	2493	Nurse Call Connection Lead	6.35mm mono to 3.5mm mono (TS) – 2m
	4151	Break cord for 4114 mat to monitor	USB A to RJ10 (wired pin 1 and 4) – 3m

cura: GUARANTEE: All Products supplied with full warranty and guaranteed to be free from defects in material and workmanship. The range of wireless devices may vary and be dependent on your specific environmental conditions. cura: products are designed to augment well designed care programs and practices, cura: does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision, and are at risk.