



Our ActiveCare LCD Pager can be standalone or work in conjunction with the ActiveCare console or 6 Channel Monitor.

It can be programmed with up to 9 devices.

Ideal for facility and home use.

3501 ActiveCare LCD Pager

Features:

- *Large easy-to-read 7-digit LCD screen displays device type and ID.
- * Can be programmed to monitor up to 9 wireless devices.
- *Up to 5 simultaneous alarms can be stored in the pager's "Call Review" memory.
- *Vibration and tone selection for hard of hearing.
- *Out of range check.
- *Wireless reset capability.
- *Direct signalling from wireless devices for use in home care environments.
- *Easily cleanable buttons for a hospital setting.
- *Typical operating range: 30 – 50m (dependent on environmental conditions).
- *Power requirements: 2 x AA batteries (included).
- *Dimensions: 83 x 55 x 26mm.
- *Housing material: Durable ABS plastic.




Range of uses:

- *The ActiveCare Pager will alert when a falls risk person is at risk or requires attention depending on the wireless devices programmed. It can be used for residential facilities or home care use.
- *The pager can be used as a standalone alert in communal areas of residential facilities for people getting up from chairs or pressing their pendants.
- *The pager can be used to extend the range of the 6 Channel Monitor (2700) or the ActiveCare Console (3500).




















Novis Healthcare
www.novis.com.au
1300 738 885

HOW TO ORDER:

Start with the 3501 ActiveCare LCD Pager and then choose up to 9 wireless programmable transmitters from the range listed below.

| Picture | Code | Description |
|-----------------------------------------------------------------------------------|------|----------------------|
|  | 3501 | ActiveCare LCD Pager |

Wireless Programmable Transmitters: Choose the devices to best suit requirements. You can choose up to 9 devices. Add devices at any stage. All pads require a 2719 or 2720 transmitter.

| Picture | Code | Description | Best Feature | Transmitter Required |
|-------------------------------------------------------------------------------------|------|----------------------------------------|----------------------------------------------|----------------------|
|  | 2630 | Wireless Chair Pad | Chair pad can be placed anywhere in the room | 2719 |
|  | 2691 | Universal Chair Pad | 18mth battery life/warranty | 2720 |
|  | 2629 | Wireless Above Mattress Bed Pad | Can detect people 40kg & under | 2719 |
|  | 3515 | Wireless Under Mattress Bed Pad | People won't even know it is there | 2719 |
|  | 2670 | Wireless PreAlert Bed Pad | Alerts when they move to the edge of the bed | - |
|  | 2690 | Universal Above Mattress Bed Pad | 18mth battery life/warranty | 2720 |
|  | 2624 | Wireless Floor Sensor Mat | Can be moved around the room as required | 2719 |
|  | 2671 | Universal Single Fold Floor Sensor Mat | 18mth battery life/warranty | 2720 |
|  | 2672 | Universal Dual Fold Floor Sensor Mat | 18mth battery life/warranty | 2720 |
|  | 2692 | Universal Floor Sensor Mat | 18mth battery life/warranty | 2720 |
|  | 2719 | Replacement Transmitter | 12-15mth battery life | - |
|  | 2720 | Replacement Universal Transmitter | 18mth battery life/warranty | - |
|  | 2721 | PIR Motion Sensor | Ideal as an out of/into room sensor | - |
|  | 3503 | Door Reed Switch | Ideal as an out of/into room sensor | - |
|  | 2588 | Pendant with Lanyard | Economical & lightweight | - |
|  | 2595 | Pendant on Wristband | Low battery warning LED | - |
|  | 2599 | SOS Button | Large button, easy to see and use | - |
|  | 3502 | Handset with Reset | Easy to use wireless handset | - |
|  | 2558 | Waterproof Call Point with Reset | Pull string & push button activation | - |