

## **molift**<sup>®</sup>

# **Molift Air Tilt**

Ceiling hoist with electric tiltable sling bar

## Molift Air Tilt

Molift Air Tilt provides a better working environment for carers thanks to the electrically tiltable sling bar. It reduces the burden on the carer and provides increased comfort for the user, as no physical effort is required and manual handling is minimized.

### Repositioning without physical effort

Molift Air Tilt allows single carers to handle demanding everyday transfers and repositioning of users safely and comfortably. Without any strain, the carer can adjust the user's height and position, simultaneously and with precision. This allows the user to be easily placed in an optimal sitting position, for example in a wheelchair.

Thanks to the open design, both users and carers are given a comfortable experience. The patient gets a feeling of openness from the wide field of view, while the carer has more opportunity for eye-contact.

## A flexible solution

Molift Air 350 is a well-established ceiling hoist in the Molift family, also available in a modified version with a tilt function. The tilt function is operated with the hand control, which makes the Molift Air Tilt very intuitive and user-friendly to use.

Molift Air Tilt and Molift Rail System are suitable in most care facilities, such as hospitals, nursing homes and home environments. The solution ensures that each step contributes to maximum function and a comfortable and ergonomic lifting situation for carers and users. Molift Air Tilt offers a flexible and cost-effective investment for the healthcare sector.



Reachable from floor with string. Reset possible from floor IP24

**Specification Molift Air Tilt Sling bar** 

Emergency stop

Protection class

Weight

Length

Width Height

Tilting range

Actuator speed

Protection class



7.3 kg

580 mm 580 mm 500-790 mm 75 degrees 12,6 mm/second IPX4

Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution

## molift

### Features



#### Easy positioning in the wheelchair

The possibility of simultaneous adjustment of the angles means that the user can be easily placed in an optimal sitting position in the wheelchair.



#### Increased comfort for the user

The sling bar is very open giving the user plenty of space around their head. This contributes to good visibility, openness and the possibility of more eye contact between the user and carer.



#### Flexible connection of the sling

The clip sling attachment provides flexibility and a convenient transfer. The sling straps are attached to a stud on the sling bar. The leg strap has a flexible connection allowing it to be adjusted according to the user's needs.



#### Quick-release coupling

Molift's sling bars have a quick-release coupling. This makes it easy and smooth to switch between different sling bars without the need for any tools.



## Slings with clip fixings



High Easy Sling

Clip fixings, Polyester, SWL 220 kg

Small	HYASMCP
Medium	HYAMEDCP
Large	HYALRGCP
Extra large	HYAXLRGCP



#### High Easy Toileting Sling Clip fixings, Polyester, SWL 220 kg

Small	DA-ST3ASM
Medium	DA-ST3AMED
Large	DA-ST3ALRG
Extra large	DA-ST3AXLRG



Etac is a world-leading developer of innovative assistive devices and patient handling equipment. Since 1973, we have been committed to improving quality of life for the individual, family members and caregivers.

> For the intended purpose of the products and latest information, please visit www.etac.com





Etac