



## **Cordless Pads & Mats**



PART NUMBER: 2624





PART NUMBER: 2630

PART NUMBER: 2719



PART NUMBER: 2629



PART NUMBER: 3515



PART NUMBER: 2670

#### Features:

- Guaranted 1 year warranty
- ☐ Simple operation and maintenance, easy to clean
- ☐ For use with cura¹ Cordless Pad Monitors
- □ No more cords between monitor and pads
- ☐ Virtually no false alarms over the lifetime of the pad
- ☐ All bed / chair pads made from non-allergenic soft PVC
- ☐ Floor mat has low profile edges to reduce tripping hazards & allow easy wheelchair access
- ☐ Floor Mat made of heavy-duty non-slip polyurethane) with the bottom material PVC without Phthalate
- ☐ Vent holes in the neck reduce the risk of moisture ingress

#### ■ 2670 PreAlert Cordless Bed Pad Unique Features:

- ☐ Provides early warning when a patient or resident sits on the edge of the bed
- □ Built-in transmitter (2719 Transmitter not required)
- ☐ Guaranteed 15 month warranty

## **Technical Specifications**

#### Cordless Bed / Chair Pads & Mats

Weight range for reliable operation: 30 - 200 Kg Must be used with a 2719 Cordless Pad / Mat Transmitter

2630 Cordless Chair Pad: 380 x 255 mm

**2629** Cordless Bed Pad: 760 x 510 mm

**2624** Cordless Mat: 1220 x 610 mm

3515 Cordless Under Matress Bed Pad: 800 x 290 mm

#### 2670 PreAlert Cordless Bed Pad

Weight range for reliable operation: 30 - 200 Kg Size: 760 x 255 mm Multi-zone detection design with options for single side or both sides sensing Centre delay zone to help prevent false alarms

#### 2719 Cordless Pad Transmitter

50 m\* range 433 MHz operating frequency Minimum 1 year battery life Incorporates status / in-use LED indication

☐ **Please Note:** Our Universal Range of Pads & Mats can also operate cordlessly





# **Cordless Pads & Mats**

### **Accessories**

	Code	Product
	2700	Six Channel Cordless Monitor
	2544 & 2549	Multiport Expanders
cua	2702	Three Channel Cordless Monitor
S []	3500	ActiveCare Console
Crita,	3501	ActiveCare Pager
cua <sup>1</sup>	3508	LED Pager
	LRCS System (Long Range Call System)	