

Multiport Expanders And Wireless Programmable Transmitters





The Multiport Expander allows wireless monitoring of your patients / residents within their rooms with no installation.

Easy to use and economical.





Out Of Room Alerts

2544 MULTIPORT EXPANDER

Features:

Out Of Chair Alerts

- *Plugs directly into the socket of your call point.
- *There is no reset button on monitor, staff just cancel at the call point.
- *Program up to 2 pads / mats and 3 other transmitters as required.
- *Can also accommodate a hardwired pad.
- *Side socket for existing nurse call pendants.
- *Low battery and pad lost indicators.
- *3 minute pause mode.
- *Operating range: 8 to 12m
- *No false alarms if low battery/pad lost indicators ignored the multiport will escalate alert and trigger the call point.





Features:

- *Installed on a bracket away from the nurse call panel so it is suitable for congested nurse call bed head panels.
- *Changeable curly cord connection lead (included) allows different wiring configurations to suit most nurse call systems.
- *All other features as above monitor 2544.

Novis Healthcare www.novis.com.au 1300 738 885

HOW TO ORDER:

Monitors: Choose the Multiport Expander to best suit your call point.

Picture	Code	Description	What's in the Box
THE STATE OF THE S	2544	Multiport Expander	1 x 2544 Monitor,
			3 x Alkaline AAA Batteries,
			1 x Quick Start Guide
	2549	Multiport Expander Bracket Mount	1 x 2544 Monitor, 1 x 3504 Mounting Bracket,
			3 x Alkaline AAA Batteries, 1 x Quick Start Guide,
			1 x Curly Cord Lead (choose from below curly cords)
	4160	Curly Cord Nurse Call Connection Lead	Right angle 6.35mm Stereo (TS)
	4161	Curly Cord Nurse Call Connection Lead	Right angle 6.35mm Stereo (TR)
	4162	Curly Cord Nurse Call Connection Lead	Right angle 6.35mm Stereo (RS)
	4163	Curly Cord Nurse Call Connection Lead	Right angle 3.5mm Stereo (TS)
	4164	Curly Cord Nurse Call Connection Lead	Right angle 3.5mm Stereo (TR)
	4165	Curly Cord Nurse Call Connection Lead	Right angle 3.5mm Stereo (RS)

Wireless Programmable Transmitters: Choose the devices to best suit requirements. You can choose up to 2 pads / mats and 3 other devices. Add devices at any stage. All pads / mats require a 2719 or 2720 transmitter.

Picture	Code	Description	Best Feature	Transmitter Required
	2630	Wireless Chair Pad	Chair pad can be placed anywhere in the room	2719
	2691	Universal Chair Pad	18mth battery life/warranty	2720
	2629	Wireless Above Mattress Bed Pad	Can detect people 30kg & under	2719
1 1 1 1	3515	Wireless Under Mattress Bed Pad	People wont even know it is there	2719
(-)	2670	Wireless PreAlert Bed Pad	Alerts when they move to the edge of the bed	-
	2690	Universal Above Mattress Bed Pad	18mth battery life/warranty	2720
	2624	Wireless Floor Sensor Mat	Can be move around the room as required	2719
-	2671	Universal Single Fold Floor Sensor Mat	18mth battery life/warranty	2720
-	2672	Universal Dual Fold Floor Sensor Mat	18mth battery life/warranty	2720
	2692	Universal Floor Sensor Mat	18mth battery life/warranty	2720
Control of the second of the s	2719	Replacement Transmitter	12-15mth battery life	-
	2720	Replacement Universal Transmitter	18mth battery life/warranty	-
	2721	PIR Motion Sensor	Ideal as an out of/into room sensor	-
	3503	Door Reed Switch	Ideal as an out of/into room sensor	-
	2588	Pendant with Lanyard	Economical & lightweight	-
•	2595	Pendant on Wristband	Low battery warning LED	-
	2599	SOS Button	Large button, easy to see and use	-
2.5	3502	Handset with Reset	Easy to use wireless handset	-
. 0,	2558	Waterproof Call Point with Reset	Pull string & push button activation	-