



The 6 Channel Wireless Monitor not only allows wireless monitoring around your facility it has many features and benefits with no installation.

Easy to use and economical.



Reassurance Pendant Monitoring



Communal Area Monitoring



Standalone or through Nurse Call System




Room Monitoring

2700 6 CHANNEL WIRELESS MONITOR

Features:

- *Monitors up to 6 wireless devices simultaneously
- * Can also accommodate a hardwired pad.
- *Caregiver secure reset option
- *Alarm memory and auto re-arm
- *Can be connected to a nurse call system (alarm tone can be set to silent)
- *Operating range: up to 30m line of sight
- *3 selectable alarm tones
- *Easy thumbwheel adjustable alarm volume
- *All safety alerts sound even if alarm volume set to silent, including low battery, AC lost, Pad signal lost , Change pad (SafeTime) & nurse call lost.
- *3-minute pause mode.
- *Battery operated or optional external plug pack.






Range of uses:

- *Monitoring wireless devices in a patient/resident's room to notify if assistance is required or of falls risk movement.
- *As a reassurance pendant alert when connected to the nurse call system, audible alarm to let residents / patients know their call has gone through.
- *In communal areas to monitor up to 6 chair pads or pendants.
- *Mobile wheelchair 'out of chair' alerts. Will emit an audible alarm to notify when a wheelchair user gets up or falls out of their chair.
- *Monitors night-time pacing, alarm can be set for 15-minute delay. This gives the person time to go to the toilet but notifies if they don't get back into bed.

Novis Healthcare
www.novis.com.au
1300 738 885

HOW TO ORDER:

Start with the 2700 6 Channel Monitor, if you would like the alerts to go through your nurse call system also order a nurse call connection lead.

Picture	Code	Description	What's in the Box
	2700	6 Channel Wireless Monitor	1 x 2700 Monitor, 3 x Alkaline AA Batteries, 1 x Soft Protective Cover, 1 x Non-Absorbent Mounting Strap, 1 x Caregiver Key, 1 x Quick Start Guide
	2460	Nurse Call Connection Lead	6.35mm stereo to 6.35mm stereo Plug (TS) – 2m
	2461	Nurse Call Connection Lead	6.35mm mono to 6.35mm stereo (TS) – 2m
	2643	Nurse Call Connection Lead	6.35mm mono to 6.35mm mono (TS) – 2m – suits most systems
	2646	Nurse Call Connection Lead	6.35mm mono to 6.35mm mono (TS) – 3m – suits most systems
	2493	Nurse Call Connection Lead	6.35mm mono to 3.5mm mono (TS) – 2m

Wireless Programmable Transmitters: Choose the devices to best suit requirements. You can choose up to 6 devices. Add devices at any stage. All pads require a 2719 or 2720 transmitter.

Picture	Code	Description	Best Feature	Transmitter Required
	2630	Wireless Chair Pad	Chair pad can be placed anywhere in the room	2719
	2691	Universal Chair Pad	18mth battery life/warranty	2720
	2629	Wireless Above Mattress Bed Pad	Can detect people 30kg & under	2719
	3515	Wireless Under Mattress Bed Pad	People won't even know it is there	2719
	2670	Wireless PreAlert Bed Pad	Alerts when they move to the edge of the bed	-
	2690	Universal Above Mattress Bed Pad	18mth battery life/warranty	2720
	2624	Wireless Floor Sensor Mat	Can be moved around the room as required	2719
	2671	Universal Single Fold Floor Sensor Mat	18mth battery life/warranty	2720
	2672	Universal Dual Fold Floor Sensor Mat	18mth battery life/warranty	2720
	2692	Universal Floor Sensor Mat	18mth battery life/warranty	2720
	2719	Replacement Transmitter	12-15mth battery life	-
	2720	Replacement Universal Transmitter	18mth battery life/warranty	-
	2721	PIR Motion Sensor	Ideal as an out of/into room sensor	-
	3503	Door Reed Switch	Ideal as an out of/into room sensor	-
	2588	Pendant with Lanyard	Economical & lightweight	-
	2595	Pendant on Wristband	Low battery warning LED	-
	2599	SOS Button	Large button, easy to see and use	-
	3502	Handset with Reset	Easy to use wireless handset	-
	2558	Waterproof Call Point with Reset	Pull string & push button activation	-