





FallsPro is a standalone system enabling the Aged Care / Hospital sector to monitor many falls risk residents / patients at their most vulnerable over a large area facility.

It is cost effective as you do not require a monitor for every room.

FallsPro System



Features:

- *Ideal for dementia wards, surgical and medical wards.
- *Monitors separate area / wards.
- *Fully wireless no cabling required for easy installation, configuration, operation, expansions, and upgrades.
- *Monitor up to 96 wireless bed and chair pads.
- *Individually identifies each resident bed or chair pad.
- *2-line backlit LCD display.
- *High visibility flashing LED & adjustable audible alarm.
- *Displays low battery and fault conditions.
- *Serial output for printer or RF paging.

Novis Healthcare www.novis.com.au 1300 738 885

REQUEST A QUOTE:

Please see below for information regarding the FallsPro monitor and its accessories.

Contact us and we will work with you to ascertain your requirements and provide you with a no-obligation quote including installation.

FallsPro System Equipment:

Picture	Code	Description	Best Features
	2501	LRCS Care Station	Up to 96 wireless bed and chair pads. Monitors up to 10 separate area ID's. 2-line backlit LCD display. High visibility flashing LED & adjustable audible alarm. Displays low battery & fault conditions. Serial output for printer or RF paging.
cura'	2509	LRCS Database Software	Connects to 2501 Care Station. Displays alarm call information. Resident/patient database. Displays low battery alerts for transmitters. Alarm & system real time log & history. System diagnostics and fault warnings. Please note: Does not include PC operation system or hardware.
	2502	Zone Repeater	Monitors up to 64 wireless call points/transmitters. Communicates directly with the 2501 Care Station. Up to 48 Zone Repeaters can work with station. LED status indications. Each zone repeater increase range by approximately 250m. Fully supervised wireless operation.
	2543	Paging Interface	Connects directly to 2501 Care Station's serial port. Connects directly to 2547 Paging Transmitter. Connects directly to serial printer for logging. LED power & status indicators.
	2547	POCSAG VHF Paging Transmitter	Connects directly to 2543 Paging Interface serial port. Sends alphanumeric alarm messages to 2548 Pagers. Synthesized crystal & 136 – 174Mhz. Output power programmable (50W to 4W). 512 or 1200 Baud selectable. 4 Hardwired alarm inputs.
	2548	POCSAG VHF Alphanumeric Pager	Receives alphanumeric alarm messages from 2547 Paging Transmitter. 4 line back lit LCD display. Memory capacity 32000 characters. Synthesized crystal & cap code programmable. Time & date message stamping.
12:00	2508	Rechargeable Battery Pack for Zone Repeaters	6V Internal lead acid battery pack to power Zone Repeater if mains power fails. Please note: includes required connection PCB assembly.
Alarm Rolom 81.	2512	LED POCSAG Display Single Sided	Displays LRCS alarms. 1 line 10 characters scrolling display. 16 Message memory. Multicoloured LEDS. Multiple capcodes.
filarin Rolom 81	2513	LED POCSAG Display Double Sided	Displays LRCS alarms. 1 line 10 characters scrolling display. 16 Message memory. Multicoloured LEDS. Multiple capcodes.

FallsPro Wireless Bed And Chair Pads:

Every bed or chair pad requires a 2719 transmitter.

Picture	Code	Description	Best Features
	2719		Communicates directly with the 2502 Zone Repeater.
		Cordless Pad	Works with cordless bed /chair pads.
		Transmitter	Range up to 50m.
			1 year battery life.
	3515		Dimensions: 800 x 290mm. Works underneath mattress. Easy to clean non-allergenic soft PVC.
		Cordless Under	Vent holes in neck reduce the risk of moisture ingress. Recommended weight range: 40-200kg.
		Mattress Bed Pad	Requires a 2719 transmitter.
			Guaranteed 12-month warranty.
	2629		Dimensions: 760 x 510mm. Easy to clean non-allergenic soft PVC.
		Cordless Above	Vent holes in neck reduce the risk of moisture ingress. Recommended weight range: 30-200kg.
		Mattress Bed Pad	Requires a 2719 transmitter.
			Guaranteed 12-month warranty.
	2630		Dimensions: 380 x 255. Easy to clean non-allergenic soft PVC.
		Cordless Chair	Vent holes in neck reduce the risk of moisture ingress. Recommended weight range: 30-200kg.
		Pad	Requires a 2719 transmitter.
			Guaranteed 12-month warranty.