



Newport Floorline Bed

A fully-featured electric bed designed for comfort, safety and versatility

Features


Integrated Bed Extension

Built-in length and width extension options provide a 100-200 mm adjustment. The extension kit includes a wipe-clean mattress platform that smoothly slides out at the foot end, offering additional space and flexibility to accommodate varying patient needs.



One Touch Chair Positioning

The auto contour feature is designed to simultaneously raise both the head and knees, ensuring optimal positioning.




Safety Descent

The bed's lowering speed slows prior to reaching its lowest height, allowing carers additional time to remove obstacles as the bed gently lowers to the floor.


Auto-Regression

The backrest auto-regression mechanism swings up and back as the bed profiles for comfort and reduces shear and friction forces that can assist in the prevention of skin damage.



Underbed Lighting

The Underbed Light is a standard integrated safety feature designed for enhanced visibility and security.



Linked Braking

Easily lockable castors at the foot end allow for secure, single-motion locking. This system ensures the bed remains stable and stationary during use. When raising or lowering, the locking mechanism prevents unwanted movement, eliminating the risk of footprint creep or damage to nearby walls.

Backlit Hand Controller

Hand controller features clearly labeled icons and individually lockable buttons, ensuring ease of use and safety.

Accessories



Back-Up Battery

BE9013

A rechargeable battery pack designed to ensure reliable operation during power failures. The battery conveniently attaches beneath the bed frame, providing continuous functionality when needed.



Self-Help Pole

BE9014

A support pole designed to assist with patient repositioning while lying in bed. It easily mounts into inbuilt brackets at either side of the bed's head.



Attendant Control Panel

BE9030

Allows caregivers functional control and the ability to lock out certain bed functions from being operated via the handset.

Clamp On Support Bed Rails

The clamp on support bed rails are designed for maximum patient and carer safety, featuring a tool free attachment system that provides extra assistance to help the user get in and out of bed. [Available in 3 sizes.](#)


Clamp On Grab Bar (400 mm)

BE9035




Long Clamp On Rail

BE9036



Narrow Clamp On Rail

BE9037



3/4 Rail

BE9017 (sold in pairs)

3/4 length tubular steel bed rails are suitable for use on Newport bed. The bed rails bolt securely on to the bed frame and can be folded down out of the way.



Specifications

	Single	King Single
Code	BE1010	BE2010
Safe working load	225 kg	250 kg
Capacity	190 kg	215 kg
Warranty	10 year bed frame 5 year electrical 1 year castor	
Frame material	Powder coated steel	
Head / foot panel	Colour	Oak Laminate (standard)
	Style	Contemporary
Platform height range	110 to 740 mm	
Extension length	200 mm	
Bed frame weight	115 kg	120 kg
Dimensions	Length	2115 to 2310 mm
	Width	920 mm
Mattress size recommended	Length	1980 to 2000 mm
	Width	880 to 900 mm
	Height	150 mm
Backrest angle	0 to 70°	
Knee break angle	0 to 35° (4 Section Bed)	
Trendelenburg tilt	0 to 15°	
Compliance	AS/NZS IEC 60601-2-52	

Accessories

BE9013	Backup Battery
BE9014	Self Help Pole
BE9015	IV Pole
BE9017	3/4 Side Rail, Pair
BE9030	Attendant Control Panel
BE9035	Clamp on Grab Bar (400 mm)
BE9036	Long Clamp on Rail
BE3037	Narrow Clamp on Rail



Find out more about the Newport Bed at
www.novis.com.au/product/newport-floorline-bed

For custom headboard or footboard configurations, please reach out to your Novis representative for assistance.